TERMOGRAVIMETRIK ANALİZÖR

TGA THERMOSTEP

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FAYDALARI

- | Measurement of up to 19 samples in one analysis
- | Sample weights of up to 5 g
- | Fast heating rates, accurate temperature control
- | High-performance, precise weighing cell
- | Automatic placing and lifting of crucible covers
- | Robust design allows for use in laboratories and production


TERMOGRAVIMETRIK ANALIZÖR THERMOSTEP TGA

RELIABLE AND FLEXIBLE

ELTRA TGA analyzers are an ideal alternative to standard laboratory ovens and muffle furnaces for thermogravimetric analysis. Thanks to a programmable furnace that is connected to an integrated balance, heating and weighing are combined in one instrument. This saves time-consuming manual work and allows for high sample throughput. In addition, typical parameters such as moisture, ash and volatiles can be determined in one analysis run.

The TGA Thermostep processes up to 19 different samples, typically weighing between 500 mg and 5 g, in one analysis cycle. The surrounding atmosphere and temperature of up to 1,000 °C within the heating chamber can be freely defined by the user during analysis to create a standard operating procedure. The crucible lids, covering of the samples, can be raised or lowered at each stage of the analysis, thus allowing for safe and ASTM-compliant determination of volatiles in coal samples.

Encapsulated weighing cell

PRECISE RESULTS

HIGH-PERFORMANCE ANALYSIS TECHNOLOGY

The TGA Thermostep is a powerful thermogravimetric analyzer characterized by robust design, high precision and flexibility. It is possible to apply different atmospheres and to use sample weights of up to 5 g. The Thermostep reliably and efficiently measures parameters such as moisture, ash and volatiles according to a user-defined SOP.

PURGING GAS

The TGA Thermostep is very flexible with regards to the purging gas used. At each stage of the analysis either nitrogen, oxygen or surrounding atmosphere can be selected. In the latter, the surrounding atmosphere penetrates into the TGA Thermostep, gently oxidizing the samples.

TEMPERATURE CONTROL

The furnace temperature is monitored by two thermocouples which are not encapsulated. One thermocouple monitors the temperature inside the furnace, the other monitors the temperature within the heating element. Due to the absence of the encapsulation the heating can be controlled quickly and precisely.

ENCAPSULATED WEIGHING CELL

The latest TGA Thermostep generation features an encapsulated weighing cell with 0.1 mg resolution providing highly precise measurements. The encapsulation isolates the weighing cell from the ambient atmosphere and is extremely stable. The weighing cell is connected to the furnace by a ceramic spindle with pedestal on which the crucibles are placed.

HIGH-CAPACITY HEATING ELEMENTS

COOLING

TERMOGRAVIMETRIK ANALIZÖR THERMOSTEP TGA

SIMPLE OPERATION YIELDS QUICK RESULTS

Operation of the TGA Thermostep is simple, convenient and safe. After selecting the Standard Operating Procedure (SOP) in the PC, the sample ID's can be entered into the software. The samples are then weighed in the crucible at the position assigned to the sample ID in the carousel. After one sample has been weighed, the carousel automatically rotates to the next position and the next registered sample can then be weighed in the crucible. Alternatively, a carousel filled with samples which has been weighed externally, can be placed into the analyzer.


It is also possible to position a second carousel with crucible lids above the crucibles. Once the analysis is finished, a new cycle can be started after a short cool-down period.

Weighing the sample

Option: Crucible lids


Display of analysis results

TERMOGRAVIMETRIK ANALİZÖR THERMOSTEP TGA

PC CONTROL WITH WINDOWS[®]-BASED SOFTWARE

The ELTRA TGA Thermostep software provides all necessary features for a fast, easy, safe and reliable thermogravimetric analysis.

ANALYSIS GRAPH

The TGA software can provide different content of the TGA measurement

- | Weight loss graph
- | Absolute weights
- | First deviation of weight loss
- | Temperature graph
- | Measuring stages


ANALYSIS RESULTS

Time	Weight	Temp	Rate	dt	dt%
0:00	1.000	1000	0.000	0.000	0.000
0:05	0.995	1000	-0.005	0.000	-0.500
0:10	0.980	1000	-0.015	0.000	-1.500
0:15	0.950	1000	-0.030	0.000	-3.000
0:20	0.900	1000	-0.050	0.000	-5.000
0:25	0.850	1000	-0.050	0.000	-5.000
0:30	0.800	1000	-0.050	0.000	-5.000
0:35	0.750	1000	-0.050	0.000	-5.000
0:40	0.700	1000	-0.050	0.000	-5.000
0:45	0.650	1000	-0.050	0.000	-5.000
0:50	0.600	1000	-0.050	0.000	-5.000
0:55	0.550	1000	-0.050	0.000	-5.000
1:00	0.500	1000	-0.050	0.000	-5.000
1:05	0.450	1000	-0.050	0.000	-5.000
1:10	0.400	1000	-0.050	0.000	-5.000
1:15	0.350	1000	-0.050	0.000	-5.000
1:20	0.300	1000	-0.050	0.000	-5.000
1:25	0.250	1000	-0.050	0.000	-5.000
1:30	0.200	1000	-0.050	0.000	-5.000
1:35	0.150	1000	-0.050	0.000	-5.000
1:40	0.100	1000	-0.050	0.000	-5.000
1:45	0.050	1000	-0.050	0.000	-5.000
1:50	0.000	1000	-0.050	0.000	-5.000
1:55	0.000	1000	-0.050	0.000	-5.000
2:00	0.000	1000	-0.050	0.000	-5.000

The results table displays the current available measurement data.

- | Display of new results after every measuring step
- | Results can be calculated as analysed or on dry base
- | Optional grouping for multiple measurements of one sample type

USER DEFINED SCREEN LAYOUT

For a better overview on the current analysis status, the user can create an own layout of the screen

- | Measuring graphs
- | Current results
- | Position of the carousel
- | Diagnosis screens

TERMOGRAVİMETRİK ANALİZÖR THERMOSTEP TGA

UYGULAMALAR

Due to the flexible application management the TGA Thermostep can measure TGA parameters in many samples like coal, Kok, Yapı malzemeleri, Gıda ...

ANALYSIS OF COAL & COKE

The TGA Thermostep can determine the parameter moisture, volatile and ash in one analysis cycle due the automatic placing and removing of lids. The TGA is compliant to ASTM D 7582.

ANALYSIS OF BUILDING MATERIALS

Building materials require the measurement of moisture and loss on ignition (LOI). The LOI can indicate the presence of organic matter or carbonates.

ANALYSIS OF FOOD

Some foods like pet foods, biscuits, or milk require the measurement of dry mass and ash content. Due to the encapsuled weighing cell even the ash content of flour could be measured safe and precise.

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